## **SMALL BATCH** LABORATORY GRINDING



Simple, economical and capable of producing small lot dispersions



L-1 through L-3-J Models can be adapted to a drill press with a minimum 1/2hp motor and a spindle speed of 2,400 RPM



Optional power stand for the L-1-through L-3-J with drive motor and starter package in either explosion-proof Class I, Div 1, Group D or non-explosion-proof configuration



All lab grinding mills have stainless steel contact parts



Agitator discs available in either cast allow or nylon



An optional screen cover is available for extended dwell times on L-2, L-3, and L-3-J Models





L-1 (1 LITER)

Water-jacketed for 300ml samples

• Allows you to determine end point on color/particle size

DIVISION OF CHICAGO BOILER COMPANY - SINCE 1891 -

• Can be used with media < 0.1mm

## L-2 (1 LITER)

Top-feed continuous mill

- 45-mesh screen for use with media as small as 0.6mm
- Allows for quick bracketing of vehicle/ pigment loading

## L-3 (2 LITER)

Same as the L-2 except for increased capacity

## L-3-J (2 LITER)

Same as L-3 except includes a water jacket for applications where temperature control is necessary



\*Optional Stainless Steel Configuration shown above